Description

MultiCore® Flow is a dual-curing, fluoride-containing, radiopaque composite for core build-ups. It cures chemically without the use of light. Light-curing is optional. MultiCore Flow is available in a double-barrel cartridge or a double-push syringe with automixing tips for fabricating core build-ups using matrices.

Composition

The monomer matrix consists of dimethacrylate (29 wt %).

The inorganic fillers are barium glass, ytterbiumtrifluoride, Ba-Al-fluorosilicate glass and highly dispersed silicon dioxide (70 wt %).

Additional contents are catalysts, stabilizers and pigments (1 wt %).

The total content of inorganic fillers is 46 vol %. The particle size ranges from 0.04 to 25 µm.

Shades

MultiCore is available in four shades: The dentin shades Light (A1/B1), Medium (A2/A3) and White are particularly suitable for use under all-ceramic restorations, and Blue is recommended for use under metal or porcelain-fused-to-metal crowns.

Working time

90 to 120 s (at 37°C/98.6°F)

Mixing ratio

The base and catalyst pastes of MultiCore Flow are mixed at a ratio of 1:1. MultiCore Flow is mixed by pressing the pastes through the static mixing tip.

Indications

- Core build-up of vital and non-vital teeth
- Adhesive cementation of glass fibre-reinforced endodontic posts

Contraindications

The use of MultiCore Flow is contraindicated

- if the operating field cannot be kept sufficiently dry, or if the stipulated application technique cannot be used;
- if the patient is known to be allergic to any of the ingredients in MultiCore Flow.

Side effects

- In individual cases, components of MultiCore Flow may lead to sensitization. MultiCore Flow should not be used in such cases.
- In vital teeth, dentinal areas close to the pulp should be covered with a suitable liner.

Interactions

Substances containing eugenol/oil of cloves may inhibit the polymerization of composite materials. Consequently, application of such materials together with MultiCore Flow must be avoided.

MultiCore® Flow

EN Instructions for Use

- Two-component core build-up material

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TR Kullanma Talimatı

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I. Fabrication of a reconstructive build-up

Apply a dentin bonding agent to the clean dentin surfaces. Ivoclar Vivadent recommends using either Adhese® Universal, Syntac® or Excite® F as the bonding agent when MultiCore Flow is used in the light-cure mode. If no light is used, the application of a dual-curing adhesive, e.g. Excite F DSC, is recommended.

These instructions describe the procedure with Adhese Universal. Please refer to the instructions of the respective bonding agents for further information.

1. Isolation

A dry operating field must be ensured. Ideally, a rubber dam should be used (e.g. OptraDam® Plus).

2. Pulp protection / Base

In very deep cavities, areas close to the pulp should be selectively coated with a calcium hydroxide liner (e.g. ApexCal®) and subsequently covered with pressure-resistant cement (glass ionomer cement, e.g. Vivaglass® Liner).

3. Conditioning with phosphoric acid gel (optional)

Apply phosphoric acid gel (e.g. Total Etch) to the prepared enamel and then flow the etchant onto the prepared dentin. The etchant should be left to react on the enamel for 15–30 s and on the dentin for 10–15 s. Then rinse thoroughly with a vigorous stream of water for at least 5 s and dry with compressed air until the etched enamel surfaces appear chalky white.

4. Application of the adhesive

- Starting with the enamel, thoroughly coat the tooth surfaces to be treated with Adhese Universal.
- The adhesive must be scrubbed into the tooth surface for at least 20 s. This time must not be shortened. Applying the adhesive on the tooth surface without scrubbing is inadequate.
- Disperse Adhese Universal with oil- and moisture-free compressed air until a glossy, immobile film layer results.
- Light-cure Adhese Universal for 10 s using a light intensity of $\geq 500 \text{ mW/cm}^2$.
- Contamination by blood or saliva must be avoided after the application of the bonding agent.

Fig. 1

Fig. 2

Fig. 3

Fig. 4

5. Insertion of the cartridge

(If you are using an automix syringe, please proceed to point c)

- Press the black release lever located below the plunger on the back of the dispenser and pull plunger back as far as possible (Fig. 1).
- Lift the cartridge lock and insert the cartridge with the "V" shape on the cartridge base turned downwards. Lower the cartridge lock (Fig. 2).
- Remove the sealing cap by turning it 1/4 turn counter clockwise.
- It is **very important** to clear or bleed the cartridge or syringe **prior** to applying the mixing tip. Gently press the dispenser lever until both components (base and catalyst) begin to flow out evenly onto a mixing pad.
- For each application, place a new mixing tip on the cartridge/syringe. Push the tip completely down until the notch of the mixing tip is aligned with the notch on the housing of the cartridge/syringe. Secure the mixing tip in place by gripping the coloured base, and not the mixing tip, and by turning it 1/4 turn clockwise (Fig. 3). As MultiCore Flow will cure in the used mixing tip, it may serve as a seal for the contents of the cartridge/syringe until needed once again (replace with a new tip just before the next use).
- If necessary apply an intra-oral tip to the mixing tip (Fig. 4). Larger amounts of material can be extruded if the intra-oral tip is shortened by approx. 1 mm.

6. Application of MultiCore Flow

- Place matrix band around the prepared tooth.
- Apply MultiCore Flow directly to the cavity with the intra-oral tip in place.
- If no light is used, wait for 4–5 minutes until the chemical curing process is complete before removing the matrix. If the material is light-cured, it may be ground immediately after completing the curing cycle.
- Optionally, MultiCore Flow may be cured from the occlusal for approx. 10 s (approx. $1,100 \text{ mW/cm}^2$, e.g. Bluephase Style) ensuring a minimal distance between the light emission window and the build-up. Remove the matrix band. If necessary, cure incompletely polymerized areas for another 10 s.
- Use rotary instruments to complete the crown preparation.

II. Adhesive cementation of endodontic posts made of fibre-reinforced composite and core build-ups

1. Pretreatment of the root canal and the coronal dentin

- A dry operating field must be ensured. Ideally, a rubber dam should be used, e.g. OptraDam Plus.
- Prior to cementing endodontic posts, clean the root canal thoroughly to remove any residue of root canal sealers (Residue of eugenol-based sealers may inhibit the polymerization of the luting composite).
- Condition the root canal, the coronal dentin and the enamel surfaces according to the instructions of the manufacturer of the dentin bonding agent used (eg or Excite F DSC).
- After applying the dentin bonding agent, make sure to remove all excess material from the root canal using paper points.
- Do not light-cure the bonding agent.
- Place matrix band around the tooth.

2. Adhesive cementation

- Clean and condition post surface according to the manufacturer's instructions.
- Coat the endodontic post, which has been prepared according to the instructions of the post manufacturer, with the mixed MultiCore Flow. Additionally, MultiCore Flow may be applied to the root canal by means of the Intra Canal Tips. It is not recommended to use a lentulo

spiral for loading MultiCore Flow into the root canal, as this will lead to the intermixing of bonding agent excess and MultiCore Flow. As a result, the working time of MultiCore Flow will be reduced, so that it may be impossible to push the post into the desired position within the root canal.

- Excess cement will be displaced when the post is inserted.
- Distribute excess cement across the occlusal preparation surface for total coverage.
- Subsequently, light-cure MultiCore Flow for 20 s (at a light intensity of approx. 1,100 mW/cm², e.g. Bluephase Style). In the process, hold the post in position with the curing light. If non-translucent, opaque endodontic posts are cemented, wait for the self-curing mechanism to take effect.

3. Core build-up

- Apply and cure MultiCore Flow as described above (see "I. Fabrication of a reconstructive build-up").

Additional notes

MultiCore Flow should be at room temperature when processed. At lower temperatures (e.g. immediately after removal from the refrigerator) it may be difficult to extrude the material.

Warnings

- Unpolymerized MultiCore Flow should not come in contact with skin, mucous membrane, or eyes. Unpolymerized MultiCore Flow may have a slight irritating effect and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

Storage

- Storage temperature: 2–8 °C/36–46 °F
Storage stability can only be guaranteed if the material is stored at this temperature.
- Leave the used mixing tip on the MultiCore Flow cartridge/syringe. It serves as a seal.
- Close syringes and cartridges immediately after use.
- Shelf life: see date of expiration on label.

Keep out of the reach of children!

For use in dentistry only!

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174, Australia

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www.ivoclarvivadent.us

Ivoclar Vivadent Shanghai

Trading Co., Ltd., 2/F Building 1, 881 Wuding Road, Jing An District

200040 Shanghai, China

Tel. +86 21 6032 1657, Fax +86 21 6176 0968

www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520, Bogotá, Colombia
Tel. +57 1 627 3399, Fax +57 1 633 1663

www.ivoclarvivadent.co

Ivoclar Vivadent SAS

B.P. 118, F-74410 Saint-Jorioz, France
Tel. +33 4 50 88 64 00, Fax +33 4 50 68 91 52

www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2, D-73479 Ellwangen, Jagst, Germany
Tel. +49 7961 889 0, Fax +49 7961 6326

www.ivoclarvivadent.de

Wieland Dental + Technik

GmbH & Co. KG, Schwenninger Strasse 13, D-75179 Pforzheim, Germany
Tel. +49 7231 3705 0, Fax +49 7231 3579 59

www.wieland-dental.com

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza, 15 B Shah Industrial Estate
Veera Desai Road, A, ndheri (West), Mumbai, 400 053, India

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www.ivoclarvivadent.com

Ivoclar Vivadent Ltd.

12 Omega St, Rosedale, PO Box 303011 North Harbour
Auckland 0751, New Zealand

Tel. +64 9 914 9999, Fax +64 9 914 9990

www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawla II 78, 00-175 Warszawa, Poland

Tel. +48 22 635 5496, Fax +48 22 635 5469

www.ivoclarvivadent.pl

Ivoclar Vivadent Marketing Ltd.

Prospekt Andropova 18 korp. 6/, office 10-06, 115432 Moscow, Russia

Tel. +7 499 418 0300, Fax +7 499 418 0310

www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St., Siricon Building No.14, 2nd Floor, Office No. 204

P.O. Box 300146, Riyadh 11372, Saudi Arabia

Tel. +966 11 293 8345, Fax +966 11 293 8344

www.ivoclarvivadent.com

Ivoclar Vivadent Pte. Ltd.

171 Chin Swee Road, #02-01 San Centre, Singapore 169877

Tel. +65 6535 6775, Fax +65 6535 4991

www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

C/ Ribera del Loira nº 46, 5ª planta, 28042 Madrid, Spain

Tel. + 34 913 757 820, Fax + 34 913 757 838

www.ivoclarvivadent.es

Ivoclar Vivadent AB

Dalvägen 14, S-169 56 Solna, Sweden

Tel. +46 8 514 939 30, Fax +46 8 514 939 40

www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

: Tesvikiye Mahallesi, Sakayik Sokak, Nisantas' Plaza No:38/2

Kat:5 Daire:24, 34021 Sisli – Istanbul, Turkey

Tel. +90 212 343 0802, Fax +90 212 343 0842

www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Ground Floor Compass Building, Feldspar Close, Warrens Business Park,

Enderby, Leicester LE19 4SE, United Kingdom

Tel. +44 116 284 7880, Fax +44 116 284 7881

www.ivoclarvivadent.co.uk

Ivoclar Vivadent, Inc.

175 Pineview Drive, Amherst, N.Y. 14228, USA

Tel. +1 800 533 6825, Fax +1 716 691 2285

www.ivoclarvivadent.us