Apexit[®] Plus

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Índice

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- 23 Al Khatar N, Kunzelmann KH, Hickel R (1999). Dichtigkeit und Wandständigkeit verschiedener endodontischer Sealer in schmierschichtfreien Wurzelkanälen. *DZZ* 54:479-483.
- 24 Kokkas AB, Boutsoukis ACh, Vassiliadis LP, Stavrianos CK: The influence of the smear layer on dentinal tubule penetration depth by three different root canal sealers: an in vitro study. *J Endod.* 2004 Feb;30(2):100-2
- 25 Limkangwalmongkol S, Burtscher P, Abbott PV, Sandler AB, Bishop BM (1991). A comparative study of the apical leakage of four root canal sealers and laterally condensed gutta-percha. *J Endod* 17:495-499
- 26 Limkangwalmongkol S, Abbott PV, Sandler AB (1992). Apical Dye Penetration with Four Root Canal Sealers and Gutta-Percha Using Longitudinal Sectioning. *J Endod* 18:535-539.
- 27 Chailertvanitkul P, Saunders WP, Mackenzie D (1996). An assessment of microbial coronal leakage in teeth root filled with gutta-percha and three different sealers. *Int Endod J* 29:387-392.
- 28 Chailertvanitkul P, Saunders WP, MacKenzie D (1997). Coronal leakage in teeth root-filled with gutta-percha and two different sealers after long-term storage. *Endod Dent Traumatol* 13:82-87
- 29 Chailertvanitkul P, Saunders WP, MacKenzie D (1997). Coronal leakage of obturated root canals after long-term storage using a polymicrobial marker. *J Endod* 23:610-613
- 30 Khademi AA, Ravandoost Y, Tabibian A (2004). Sealing ability of five root canal sealers against coronal leakage of *Enterococcus faecalis*. *Endodontic Practice*:31-34
- 31 Timpawat S, Amornchat C, Trisuwan W (2001). Bacterial coronal leakage after obturation with three root canal sealers. *J Endod* 27:36-39
- 32 Kokkas AB, Boutsoukis ACh, Vassiliadis LP, Stavrianos CK: The influence of the smear layer on dentinal tubule penetration depth by three different root canal sealers: an in vitro study. *J Endod.* 2004 Feb;30(2):100-2
- 33 Miletic I, Anic I, Pezelj-Ribaric S, Jukic S: Leakage of five root canal sealers. *Int Endod J.* 1999 Sep;32(5):415-8.
- 34 Miletic I, Ribaric SP, Karlovic Z, Jukic S, Bosnjak A, Anic I: Apical leakage of five root canal sealers after one year of storage. *J Endod.* 2002 Jun;28(6):431-2
- 35 Di Lenarda R: Evaluation of the clinical performance of Apexit Plus comparing thermafil technique versus non stop wave technique. University of Trieste
- 36 Pascon EA: Evaluation of the clinical performance of Apexit Plus using warm guttapercha. University of Toronto
- 37 Lipski M: Studies comparing the efficacy of root canal filling with gutta-percha lateral condensation and Thermafil obturators. *Ann Acad Med Stetin.* 2000;46:317-30.
- 38 Briseno MB, Kremers L: Klinische und röntgenologische Untersuchung des auf Calciumhydroxid aufgebauten Apexit-Wurzelkanal-Zementes. *Dental Spiegel* März 1992
- 39 Kastakova A, Halackova Z: Apikal leakage of Thermafil and lateral compactation. Part III: Sealer Apexit. The importance of sealer application. *Ces. Stomat.* 101, No 6 (205-211). 2001 [abstract]
- 40 Halackova Z: Thermafil – Effect of resin and hydroxide sealer on quality of root filing. *Ces Stomat* 99, No 1(830-36) 1999 [abstract]
- 41 RTECS Database
- 42 RCC: Acute Oral Toxicity Study with Disalicylat S55 in rats. Study report 076296, 1986.
- 43 RCC: Cytotoxicity assay in vitro: evaluation of materials for medical devices (XTT-Test) with Apexit II. 2004
- 44 Geurtsen W, Leinenbach F, Krage T, Leyhausen G: Cytotoxicity of four root canal sealers in permanent 3T3 cells and primary human periodontal ligament fibroblast cultures. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod.* 1998 May;85(5):592-7.

- 45 Miletic I, Anic I, Karlovic Z, Marsan T, Pezelj-Ribaric S, Osmak M: Cytotoxic effect of four root filling materials. *Endod Dent Traumatol.* 2000 Dec;16(6):287-90.
- 46 Beltes P, Koulaouzidou E, Kotoula V, Kortsaris AH. In vitro evaluation of the cytotoxicity of calcium hydroxide-based root canal sealers. *Endod Dent Traumatol.* 1995 Oct;11(5):245-9.
- 47 Vajrabhaya L, Sithisarn P, Wilairat P, Leelaphiwat S: Comparison between Sulphorhodamine-B dye staining and 51Cr-release method in cytotoxicity assay of endodontic sealers. *J Endod.* 1997 Jun;23(6):355-7.
- 48 Vajrabhaya L, Sithisarn P: Multilayer and monolayer cell cultures in a cytotoxicity assay of root canal sealers. *Int Endod J.* 1997 Mar;30(2):141-4.
- 49 Schwarze T, Fiedler I, Leyhausen G, Geurtsen W: The cellular compatibility of five endodontic sealers during the setting period. *J Endod.* 2002 Nov;28(11):784-6.
- 50 Schwarze T, Leyhausen G, Geurtsen W: Long-term cytocompatibility of various endodontic sealers using a new root canal model. *J Endod.* 2002 Nov;28(11):749-53.
- 51 RCC: Acute Oral Toxicity Study in Rats. 1988
- 52 Salmonella typhimorium reverse mutation assay (AMES) with Apexit II. Study report 811503, RCC-CCR, 2004
- 53 Schweikl H, Schmalz G: Evaluation of the mutagenic potential of root canal sealers using the Salmonella microsome assay. *J Mater Sci: Mater Med* 2:181-185
- 54 Heil J, Reifferscheid G, Waldmann P, Leyhausen G, Geurtsen W: Genotoxicity of dental materials *Mutat Res.* 1996 Jul 5;368(3-4):181-94.
- 55 Bratel J, Jontell M, Dahlgren U, Bergenholtz G: Effects of root canal sealers on immunocompetent cells in vitro and in vivo. *Int Endod J.* 1998 May;31(3):178-88.
- 56 Leimgruber R: Pharmakologisches Sachverständigen-Gutachten über Apexit, RCC 1990
- 57 Kolokouris I, Economides N, Beltes P, Vlemmas I: In vivo comparison of the biocompatibility of two root canal sealers implanted into the subcutaneous connective tissue of rats. *J Endod.* 1998 Feb;24(2):82-5.
- 58 Silva LA, Leonardo MR, Faccioli LH, Figueiredo F: Inflammatory response to calcium hydroxide based root canal sealers. *J Endod.* 1997 Feb;23(2):86-90.
- 59 Leonardo MR, Silva LA, Utrilla LS, Assed S, Ether SS: Calcium hydroxide root canal sealers--histopathologic evaluation of apical and periapical repair after endodontic treatment. *J Endod.* 1997 Jul;23(7):428-32.
- 60 Bernath M, Szabo J: Tissue reaction initiated by different sealers. *Int Endod J.* 2003 Apr;36(4):256-61.

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Ivoclar Vivadent AG