# Systemp.<sup>®</sup> inlay Systemp.<sup>®</sup> onlay

## Bruksanvisning

## Bruksanvisning

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## Bruksanvisning

## Productinformatie

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## Instructions for Use

### Description:

Systemp.inlay and Systemp.onlay are light-curing single-component materials for temporary restorations, generally not requiring the application of a temporary cement. Both materials are available in a transparent and universal shade and are supplied in syringes and Cavifils.

### Composition:

#### Systemp.inlay

The monomer matrix consists of high-molecular dimethacrylates. The fillers are highly dispersed silicon dioxide and copolymers. Triclosan, catalysts and stabilizers are additional ingredients.

#### System.onlay

The monomer matrix consists of high-molecular dimethacrylates. The fillers are highly dispersed silicon dioxide and copolymers. Triclosan, catalysts and stabilizers are additional ingredients.

### Indication:

- Systemp.inlay is especially suitable for deep inlay preparations with parallel walls. Systemp.inlay can also be used if small undercuts are present.
- Systemp.inlay is also suitable for relining prefabricated, temporary crowns and bridges made of polycarbonate or methacrylates.
- Systemp.inlay is also suitable for sealing implant screw access openings.
- Systemp.onlay is especially suitable for large preparations (onlays). If undercuts are present, please note that after curing, Systemp.onlay is less elastic than Systemp.inlay, which may complicate removal.

### Contraindication:

Unsuitable as temporary crown and bridge material. Systemp.inlay or Systemp.onlay should not remain in the oral cavity for more than one month. The product should not be used in patients with a known allergy to any of its ingredients.

### Side effects:

After long periods of service, the material may demonstrate discoloration. However, this does not compromise the performance of the material. In rare cases, Systemp.inlay or Systemp.onlay may cause a sensitizing reaction in patients with a hypersensitivity to any of the materials' ingredients. In these cases, the materials should no longer be used.

### Interactions:

Phenolic substances (e.g. eugenol) may inhibit the polymerization of Systemp.inlay or Systemp.onlay. Therefore, Systemp.inlay or Systemp.onlay must not be applied in combination with zinc oxide-eugenol cements.

### Application method:

#### Preparation appointment

1. Prepare, place base or liner and take impression as usual.
2. Dry prepared teeth; use a matrix if necessary.
3. Clean cavity with water spray.
4. Application of Systemp.desensitizer: Apply Systemp.desensitizer and brush the material gently into the dentin for 10 seconds.
5. Carefully dry with blown air.
6. Place a sufficient quantity of Systemp.inlay or Systemp.onlay into the cavity, using a spatula or another suitable instrument. If you use Cavifils, the material can be applied directly into the cavity by means of the Cavifil Injector. Contour the restoration, leaving minimum excess.
7. Increments of up to 4 mm can be cured with a standard curing light (e.g. bluephase® C5) in 10 seconds, if the light intensity is higher than 500 mW/cm<sup>2</sup>. The same depth of cure is achieved with LED lights that produce a higher light intensity (e.g. bluephase®) after an exposure time of 10 seconds in the low power mode (LOP).
8. Remove matrix.
9. Check occlusion.

#### Insertion appointment

1. Insert a suitable instrument (e.g. probe/scaler) into the temporary restoration. Remove Systemp.inlay or Systemp.onlay from the cavity.

2. Subsequently, clean the cavity (e.g. with rubber cups and cleaning paste).
3. Try in and seat the restoration.

#### Note 1:

Systemp.inlay/Systemp.onlay remains well in place in standard inlay preparations. If mechanical retention is minimal, the temporary restoration is cemented with a eugenol-free cement (e.g. Systemp®.link). If Systemp.inlay or Systemp.onlay is applied to undercuts of the proximal region, retention can be improved. In this case, do not place a matrix, but a wooden interdental wedge before applying Systemp.inlay/ Systemp.onlay.

#### Note 2:

Systemp.inlay/Systemp.onlay bonds with light-curing base/liner materials because of the similar chemical compositions. If such materials (light-curing glass ionomer cements, composites) are utilized, isolate the base with glycerine gel (Liquid Strip). This prevents the base/lining from being accidentally removed from the cavity together with Systemp.inlay/Systemp.onlay.

#### Note 3:

Use silicone finishers (Politip F) or tungsten carbide finishers for grinding and excess removal. Polish with silicone rubber polishers (Politip P). Excess material may also be removed with a scalpel.

#### Note 4:

To increase the stability of large cavities, a composite (e.g. Heliomolar® or Tetric EvoCeram®) may be applied to the occlusal part of the temporary. Light- cure Systemp.inlay/ Systemp.onlay separately or together with the restorative material.

#### Note 5:

Systemp.inlay / Systemp.onlay can be contoured more easily by wetting the instrument with Systemp.desensitizer or an unfilled bonding agent (e.g. Heliobond).

#### Warning:

Unpolymerized Systemp.inlay/Systemp.onlay should not come in contact with skin, mucous membrane, or eyes. Unpolymerized Systemp.inlay/Systemp.onlay may have a

slight irritating effect and may lead to a sensitization against methacrylates.

Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

#### Storage:

- Do not use Systemp.inlay/Systemp.onlay after the date of expiration.
- Storage temperature: 2 – 28 °C (36 – 82 °F).
- Close syringe / Cavifil immediately after use.
- Store material protected from light.
- Shelf life: see date of expiration.

**Keep out of the reach of children!  
For use in dentistry only!**

**Date information prepared: 12/2005  
Rev. 1**

#### Manufacturer:

Ivoclar Vivadent AG, FL-9494 Schaan/Liechtenstein

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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## **Ivoclar Vivadent AB**
## **Ivoclar Vivadent UK Limited**
## **Ivoclar Vivadent, Inc.**