Description

World-Cap® is a non-gamma-2, spherical / irregular-particle silver amalgam in self-activating capsules for dental restorations.

Composition

One Size-1 World-Cap Non-Gamma-2 capsule contains the following active substances:

Powder

Silver	183 mg
Tin	140 mg
Copper	132 mg

Liquid

Mercury	400 mg
---------	--------

One Size-2 World-Cap Non-Gamma-2 capsule contains the following active substances:

Powder

Silver	248 mg
Tin	189 mg
Copper	178 mg

Liquid

Mercury	550 mg
---------	--------

One Size-3 World-Cap Non-Gamma-2 capsule contains the following active substances:

Powder

Silver	355 mg
Tin	270 mg
Copper	255 mg

Liquid

Mercury	800 mg
---------	--------

Length of the capsule: 31.6 mm

Largest outer diameter: 13.6 mm

Indication

Class I and II restorations

Contraindication

World-Cap Non-Gamma-2 is contraindicated

- if proximal or occlusal contacts with other metal restorations are present
- if the patient suffers from impaired renal function
- if the patient is known to be allergic to any of the components of amalgam
- for retrograde or endodontic restorations
- as a core build-up material under crowns or inlays
- as a filling material under metal crowns
- for children under fifteen years of age
- for pregnant or nursing women
- if the stipulated working technique cannot be applied.

Side effects

- In individual cases, hypersensitivity, allergic reactions, or electrochemically-induced local reactions may occur.
- Mercury may have a sensitizing effect on the skin and lungs, or act as a nephrotoxin or neurotoxin.
- Electrochemical processes may cause lichen planus of the oral mucous membrane.
- After the placement or removal of amalgam restorations, a short-term increase of the mercury concentration in the blood or urine may occur.

Interactions

Mercury that has escaped during condensing or amalgam that has not set may cause amalgamation or galvanic reactions (electrochemically-induced local malaise, i.e. alterations in the sense of taste) when the restoration comes in contact with other metal restoratives. If these symptoms persist, the amalgam should be replaced with another material.

Application

Relative isolation with e.g. cotton rolls.

Moisture contamination of the amalgam before it has set completely can detrimentally affect physical properties such as stability and corrosion resistance.

Mixing times for World-Cap Non-Gamma-2

The mixing time in the Silamat® S6 is 10 s.

Mixers of other competitors should demonstrate a high starting torque and achieve approx. 4,500 rpm. The use of low-speed amalgamators is not recommended. Longer mixing times result in a softer consistency and slightly improved stability, without significant effect on the other properties.

The triturated amalgam should be bright and have a plastic consistency.

Base/liner

Cover very deep areas of the cavity near the pulp with a small amount of calcium hydroxide cement. The base/liner must be able to withstand condensation pressure. Line the cavity with a mechanically stable cement such as Vivaglass Liner (light-curing glass ionomer cement).

Placing the restoration

After mixing, open the capsule, apply the pliable amalgam to the cavity in portions and condense and model it. Ultrasonic condensers should not be used. The working time of World-Cap is approx. 4 minutes.

Polishing

The restoration should not be polished until at least 24 hours after placement. Suitable tungsten-carbide finishers or other polishing instruments (Politip-F and Politip-P) can be used. Make sure to provide steady, adequate spray cooling and suction.

Safety notes

- Mercury reacts with and embrittles particular metals and their alloys. Avoid unnecessary contact between mercury and those materials (and their alloys).
- Mercury presents a health hazard if handled incorrectly. Mercury is toxic by vapour inhalation and the effect is cumulative. Spillages of mercury should be removed immediately, including from places which are difficult to access. Use a plastic syringe to draw it up. Smaller quantities can be

World-Cap® Non-Gamma-2

EN Instructions for Use
- Non-gamma-2, 40.3% silver amalgam alloy with spherical / irregular particles in self-activating capsules

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covered by sulphur powder and removed. Individual small droplets can be picked up by tin foil (Sn) and removed. Avoid inhalation of the vapour. During and after clean-up thoroughly ventilate the area where the spill has occurred. Vacuum cleaners should not be used.

- If moisture is introduced into the amalgam before it has set, properties such as strength and corrosion resistance can be affected adversely. Use a dry field, whenever possible.
- Waste material and all primary containers that have held mercury shall be disposed of following appropriate management practice.
- As the removal of existing amalgam restorations results in the additional release of mercury, clinically sound amalgam restorations should not be removed. The removal of existing amalgam restorations should be avoided during pregnancy or lactation in particular.
- When removing amalgam restorations, the corresponding safety notes must be observed: Wear nose/mouth protection and goggles; use spray aspirators, saliva ejectors and provide adequate spray cooling.
- Only store amalgam and place and remove amalgam fillings in well-aired practice facilities.
- Avoid skin contact with unset material.
- Never open capsules prior to mixing!
- Single-use only

Storage

- **Store amalgam waste in containers that can be tightly sealed.**
- **Dispose of used capsules in accordance with national regulations.**
- Store packages at room temperature (2-28 °C/36-82 °F).
- See expiration date for shelf life.

Keep material out of children's reach!

For dental use only!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Ivoclar Vivadent Marketing Ltd.
Ivoclar Vivadent AB
Ivoclar Vivadent Limited
Ivoclar Vivadent, Inc.