### Instructions for Use

### Mode d'emploi

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### Kullanma Talimatı

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CE 0123

For dental use only.  
Caution: U.S. Federal law restricts  
this device to sale by or on the  
order of a licensed dentist.

**Manufacturer**  
Ivoclar Vivadent AG  
FL-9494 Schaan/Liechtenstein  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

ivoclar  
vivadent®  
clinical

64998832013-02-01/Rev. 0/6spr.

### Description

Multilink® N is a self-curing luting composite with light-curing option for the adhesive luting of indirect restorations made of metal, metal-ceramic, all-ceramic and composite. Multilink N is applied together with the self-etching and self-curing Multilink N Primer. The special filler composition gives Multilink N a very high radiopacity. Monobond N is recommended as a coupling agent to achieve a strong bond to precious and non-precious alloys, as well as to all-ceramics made of zirconium and aluminium oxide and silicate ceramics.

### Shades

Multilink N is available in two shades:

- transparent
- yellow

### Working time

The working and setting times depend on the ambient temperature. Once Multilink N has been dispensed from the automix syringe, the following times apply:

	At room temperature 23 °C (73 °F)	Intraorally (in combination with Multilink N Primer) approx. 37 °C (99 °F)
Working time	approx. 3 min	approx. 2 min
Curing time (including working time)	approx. 8 min	approx. 5 min

### Mixing ratio

Multilink N is always dispensed from the automix syringe in the optimum ratio. Multilink N Primer A and Multilink N Primer B are mixed in a 1:1 ratio (e.g. 1 drop Primer A + 1 drop Primer B).

### Composition

#### Multilink N

The monomer matrix is composed of dimethacrylate and HEMA. The inorganic fillers include barium glass, ytterbium trifluoride and spheroidal mixed oxide.

The particle size is 0.25–3.0 µm. The mean particle size measures 0.9 µm. The total volume of inorganic fillers is approximately 40 %.

#### Multilink N Primer A and B

Multilink N Primer A is an aqueous solution of initiators. Multilink N Primer B contains HEMA, phosphonic acid and methacrylate monomers.

### Indication

Multilink N and Multilink N Primer are used for the permanent cementation of indirect restorations where a strong bond is desired:

- Inlays, onlays, crowns, bridges and root posts made of
  - metal and metal-ceramics
  - all-ceramics, in particular opaque zirconium oxide ceramics
  - composites and fibre-reinforced composites

## Contraindication

The use of Multilink N and Multilink N Primer is contraindicated

- if a dry working field cannot be established or the stipulated working technique cannot be applied;
- if a patient is known to be allergic to any of the ingredients of Multilink N and Multilink N Primer A/B.

## Side effects

Systemic side effects are not known to date. In individual cases, allergic reactions to single components have been reported.

## Interactions

Phenolic substances (e.g. eugenol, wintergreen oil) inhibit polymerization.

Consequently, the application of products containing these components, e.g. mouth rinses and temporary cements, must be avoided.

Alkaline jet media impair the effect of Multilink N Primer.

Hydrogen peroxide and other disinfectants with an oxidative effect may interact with the initiator system, which in turn can impair the curing process. Therefore, the automix syringe must not be disinfected with oxidative agents. The syringe can be disinfected, e.g. by wiping with medical alcohol.

## Application

### 1. Removal of the temporary restoration and cleaning of the cavity

Remove any residue of the temporary luting cement from the cavity or preparation with a polishing brush and an oil and fluoride-free cleaning paste (e.g. fluoride-free Proxyl). Rinse with water spray. Subsequently, dry with air that is free of water and oil. Avoid over drying.

**Note:** Cleaning with alcohol could lead to dehydration of dentin.

### 2. Try-in of the restoration and drying

Next, check the shade, fit and occlusion of the restoration. Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently cemented, as there is a risk of fracture. If necessary, make adjustments with fine diamonds at medium speed and with slight pressure using copious amounts of water. Polish ground surfaces. When an adhesive luting protocol with composites is used, safe isolation of the operating field – preferably with a rubber dam, e.g. OptraDam, or alternatively with cotton rolls and a saliva ejector – is required. If the dental hard tissue is contaminated with blood or saliva during the try-in procedure, it must be cleaned again as described in Item 1.

### 3. Pre-treatment of the restoration

3.1 If a restoration that has been **pre-treated** in the dental lab comes in contact with blood or saliva during the try-in procedure, the contaminated bonding surface of the restoration must be cleaned in the following way:

- After try-in, thoroughly rinse the restoration with water spray and dry it with oil-free air.
- Completely coat the bonding surface of the restoration with Ivoclean using a microbrush or brush.

– Allow Ivoclean to react for 20 seconds. Thoroughly rinse the restoration with water spray and dry it with oil-free air.

- Use Monobond N on restorations made of lithium disilicate glass-ceramics or zirconium oxide ceramics (e.g. IPS e.max Press/CAD IPS e.max ZirCAD). Apply Monobond N to the cleaned surfaces using a brush or microbrush. Allow Monobond N to react for 60 seconds and subsequently disperse it with a strong stream of air.

3.2 Restorations that have not been pre-treated by the dental lab should be rinsed with water spray and dried after the try-in procedure. Thereafter, the restoration is conditioned as follows:

#### 3.2.1.1 Glass-ceramic restorations (e.g. IPS Empress®)

- Etch the restoration with 5% hydrofluoric acid (e.g. IPS® Ceramic etching gel) for 60 seconds or according to the instructions of the manufacturer of the restorative material.
- Thoroughly rinse the restoration with water spray and dry it with oil-free air.

#### 3.2.1.2 Lithium disilicate glass-ceramic restorations (e.g. IPS e.max® Press, IPS e.max CAD)

- Etch the restoration with 5% hydrofluoric acid (e.g. IPS Ceramic etching gel) for 20 seconds or according to the instructions of the manufacturer of the restorative material.
- Thoroughly rinse the restoration with water spray and dry it with oil-free air.

#### 3.2.1.3 Zirconium oxide (e.g. IPS e.max ZirCAD) and aluminium oxide ceramic restorations

- Sandblast the inner restoration surfaces (sandblasting parameters according to the instructions of the manufacturer of the restorative material).
- If necessary, clean the restoration in an ultrasonic unit for about 1 minute.
- Thoroughly rinse the restoration with water spray and dry it with oil-free air.
- **IMPORTANT!** In order to create a strong bond, do not clean the zirconium oxide surfaces with phosphoric acid.

#### 3.2.1.4 Metal or metal-supported restorations

- Sandblast the inner restoration surfaces (sandblasting parameters according to the instructions of the manufacturer of the restorative material) until a uniformly mat surface is achieved.
- If necessary, clean the restoration in an ultrasonic unit for about 1 minute.
- Thoroughly rinse the restoration with water spray and dry it with oil-free air.
- **IMPORTANT!** In order to create a strong bond, do not clean the metal surfaces with phosphoric acid.

#### 3.2.1.5 Composite or fibre-reinforced composite restorations

- Sandblast the inner restoration surfaces (sandblasting parameters according to the instructions of the manufacturer of the restorative material).
- If necessary, clean the restoration in an ultrasonic unit for about 1 minute.
- Thoroughly rinse the restoration with water spray and dry it with oil-free air.

3.2.2 Next, apply Monobond N with a brush or microbrush to the pre-treated surfaces. Allow the material to react for 60 seconds. Disperse it with a strong stream of air.

## 4. Mixing of Multilink N Primer A and Multilink N Primer B

The two primer liquids Multilink N Primer A and B are mixed in a 1:1 ratio (e.g. 1 drop Primer A and 1 drop Primer B). The mixed Primer A/B is solely self-curing and does not need to be protected against light. However, it must be applied within 10 minutes.

## 5. Application of the mixed Multilink N Primer A/B to the enamel and the dentin

Apply the mixed Multilink N Primer A/B onto the entire bonding surface using a microbrush, starting with the enamel surface, and scrub it in for 30 s.

Disperse excess with blown air until the mobile liquid film is no longer visible.

As the Primer is solely self-curing, no light-curing is necessary!

**Note:** Previous conditioning of enamel with phosphoric acid is required for Maryland bridges.

## 6. Application of Multilink N to the restoration

For each application, place a new automix tip on the syringe. Dispense Multilink N from the automix syringe and apply the desired amount onto the restoration.

As the luting material will cure in the used mixing tip, it serves as a seal for the remaining contents of the syringe until needed again (replace with a new tip before the next use).

**Note:** Multilink N should be used quickly after it has been dispensed and the restoration should be placed rapidly! Multilink N is applied directly to the inner surface of the restoration. The application of Multilink N directly to the preparation or the cavity, which has been pre-treated with Multilink N Primer, is not indicated as this would lead to considerable acceleration of the curing process and thus may compromise the fit of the restoration.

## 7. Placement of the restoration and removal of excess cement a) solely self-curing

Seat the restoration in place and secure it. Remove excess material immediately with a microbrush/brush/foam pellet/dental floss or scaler. Remove excess rapidly in difficult-to-reach areas (proximal areas, gingival margins, pontics). Due to the reaction between Multilink N and Multilink N Primer A/B, a high bond strength and a high degree of cure is achieved within a few minutes of seating the restoration.

## b) self-curing with additional light-curing (quarter technique, indicated for cases with up to 2 bridge abutments = 3 to 4-unit bridges)

Light-curing of all excess cement and subsequent removal: Seat the restoration in place and secure it. Light-cure excess cement in quarter segments (mesio-oral, disto-oral, mesio-buccal, disto-buccal) using a polymerization light at a distance of max. 10 mm. In the process, the following parameters must be observed:

Light intensity	Exposure time per quarter segment	
approx. 650 mW/cm <sup>2</sup>	3 s	e.g. Bluephase in LOW POWER mode
approx. 1,000 mW/cm <sup>2</sup>	1-2 s	e.g. Bluephase Style or Bluephase in HIGH POWER mode

Thereafter, excess cement is easy to remove with a scaler. Excess must be removed rapidly in difficult-to-reach areas (proximal areas, gingival margins, pontics). Subsequently, light-cure all margins again for 20 s (approx. 1,000 mW/cm<sup>2</sup>, e.g. Bluephase in HIGH POWER mode or Bluephase Style).

## c) self-curing with additional light-curing (quarter technique, indicated for cases with up to 6 bridge abutments = circular bridges)

Light-curing of all excess cement and subsequent removal: Seat the restoration in place and secure it. Successively light-cure excess cement on all the abutments in quarter segments (mesio-oral, disto-oral, mesio-buccal,

disto-buccal) using a polymerization light at a distance of max. 10 mm. In the process, the following parameters must be observed:

Light intensity	Exposure time per quarter segment	
approx. 650 mW/cm <sup>2</sup>	1 s	e.g. Bluephase in LOW POWER mode

Thereafter, excess cement is easy to remove with a scaler. Excess must be removed rapidly in difficult-to-reach areas (proximal areas, gingival margins, pontics). Subsequently, light-cure all the margins again for 20 s (approx. 1,000 mW/cm<sup>2</sup>, e.g. Bluephase in HIGH POWER mode or Bluephase Style).

## Note

Like all composites, Multilink N is subject to oxygen inhibition: In other words, the surface layer (approx. 100 µm) does not polymerize during the curing process, as it comes in contact with atmospheric oxygen. In order to avoid this, we recommend covering the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after the removal of excess. After complete polymerization, the glycerine gel/air block is rinsed off with water.

## 8. Finishing of the restoration

- Finish the cement joints with preparation diamonds if necessary.
- Polish restoration margins/cement joints with suitable polishers (e.g. Astropol®) as well as finishing and polishing strips.
- Check occlusion and function and make adjustments if necessary.

## Special instructions for the cementation of endodontic posts

1. For the cementation of endodontic posts, carefully clean the root canal to remove any residue of root canal fillers. (Residue of eugenol-based sealers may inhibit the polymerization of the luting composite). Ideally, the mixed Multilink N Primer A/B should be applied to the root canal and the occlusal surface of the preparation for about 15 s each using the thin (violet) microbrushes. Remove excess cement from the root canal using paper points.
2. Coat the endodontic post, which has been prepared according to the instructions of the post manufacturer, with the mixed Multilink N cement. **IMPORTANT!** Do not apply Multilink N cement with the lute into the root canal that has been wetted with Primer A/B. Otherwise, premature curing might occur, which could possibly prevent the post from being properly positioned.
3. Place the endodontic post. Excess cement will be displaced.
4. Distribute excess cement across the occlusal preparation surface for total coverage. Then, light-cure Multilink N for 20 s.
5. Core buildup material (e.g. MultiCore®) can then be directly applied onto Multilink N to serve as a bonding agent and cured according to the instructions of the manufacturer (light-curing).

## Warning

Unpolymerized Multilink N pastes and the Multilink N Primer liquids may cause slight irritation. Avoid contact with the skin, mucous membrane and eyes. If the material comes into contact with the eyes, immediately rinse with plenty of water and seek medical advice. If the material comes into contact with the skin, rinse with copious amounts of water. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

**Storage**

- Multilink N and Multilink N Primer must not be used after the expiry date.
- Storage temperature 2-8 °C / 36-46 °F.
- Multilink N Primer bottles must be tightly closed with the correct cap after use.
- In order to ensure a tight seal of Multilink N syringes, the mixing tip is left on the syringe after use.
- Expiry date: See information on bottles, syringes and packaging.

**Keep material out of children's reach.****For use in dentistry only.**

**Date information prepared:** 2013-02-01/Rev. 0

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions.

**The Material Safety Data Sheet is available online at**  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

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**Ivoclar Vivadent Limited**
**Ivoclar Vivadent, Inc.**