**Description**

Multilink® Speed is a self-adhesive, self-curing resin cement with light-curing option for the cementation of indirect restorations made of high-strength all-ceramics, metal-ceramics and metal. When Multilink Speed is used, additional dentin bonding agents are not required. The automix syringe enables time-saving direct application of the luting material into the restoration.

**Shades**

Multilink Speed is available in the shade transparent.

**Working and setting times**

The following times apply as soon as Multilink Speed has been dispensed from the automix syringe (the working and setting times depend on the ambient temperature):

	<b>At room temperature approx. 23°C/73°F</b>	<b>Intra-orally</b>
<b>Working time</b>	approx. 2 min	approx. 1 min
<b>Setting time (incl. working time)</b>	approx. 6 min	approx. 3 min

**Note**

After Multilink Speed has been dispensed from the automix syringe, the curing process can be accelerated by intensive illumination (operating light, ambient light).

**Composition**

The monomer matrix is composed of dimethacrylates and acidic monomers. The inorganic fillers are barium glass, ytterbium trifluoride, co-polymer and highly dispersed silicon dioxide. Additional contents are initiators, stabilizers and colour pigments (< 1 %). The primary particle size of the inorganic fillers is between 0.1 and 7 µm. The mean particle size is 5 µm. The total content of inorganic fillers is approx. 40 % vol.

**Indication**

- 1) Permanent cementation on natural teeth with adequate retentive features in conjunction with indirect restorations made of:
  - oxide ceramics, e.g. zirconium oxide (crowns, bridges, endodontic posts)
  - metal and metal-ceramics (inlays, onlays, crowns, bridges, endodontic posts)
  - lithium disilicate glass-ceramics (crowns, bridges)
  - fibre-reinforced composites (endodontic posts)
- 2) Permanent cementation of crowns and bridges on implant abutments that are made of the following materials:
  - oxide ceramic (e.g. zirconium oxide)
  - metal (e.g. titanium)
  - lithium disilicate glass-ceramic

**Contraindication**

Multilink Speed is contraindicated

- in situations where the preparation/shape of the implant abutment does not provide adequate retention (e.g. veneers, short or severely tapered tooth preparations/implant abutments);
- if a patient is known to be allergic to any of the ingredients of Multilink Speed.
- In general, Multilink Speed should not be applied to the exposed pulp or to dentin which is close to the pulp.

**Side effects**

Systemic side effects are not known to date. In rare cases, allergic reactions to single components have been reported.

**Interactions**

Phenolic substances (e.g. eugenol, wintergreen oil) inhibit polymerization. Consequently, the use of products, e.g. mouth rinses and temporary cements, which contain these components should be avoided. Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process. Therefore, the preparation and syringe must not be disinfected using oxidative agents. The syringe can be wiped e.g. with customary disinfecting wipes. Alkaline jet media may compromise the effect of Multilink Speed.

**Usage**

Please also refer to the relevant instructions for use of the products used in conjunction with Multilink Speed for more detailed information.

**1 Pre-treatment of the prepared tooth or implant abutment****1.1 Pre-treatment of the prepared tooth****Removal of the temporary restoration and thorough cleaning of the prepared tooth**

Remove possible residue of the temporary luting composite from the prepared tooth with a polishing brush and an oil- and fluoride-free cleaning paste (e.g. Proxyt® fluoride-free). Rinse with water spray. Then lightly dry with water- and oil-free air. Avoid overdrying.

**Note:** Cleaning with alcohol can lead to overdrying of the dentin.

**Try-in of the restoration and isolation**

Next, the shade, accuracy of fit and occlusion of the restoration can be checked.

Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently seated, because there is a risk of fracture. It is essential to ensure reliable isolation of the operating field with e.g. OptraGate®, cotton rolls, saliva ejector and absorbent pads. Dental hard tissue that has become contaminated with saliva needs to be cleaned again.

**1.2 Pre-treatment of the implant abutment****Removal of the temporary restoration**

Remove the temporary restoration and, if present, the temporary abutment. Rinse the implant lumen and peri-implant gingiva.

**Insertion and inspection of the final abutment**

According to the manufacturer's directions.

A retraction cord may be placed to be better able to check the accuracy of fit of the restoration and remove excess luting composite.

# Multilink® Speed

<b>EN Instructions for Use</b>	Self-adhesive self-curing dental resin cement with light-curing option
<b>FR Mode d'emploi</b>	Composite de collage auto-adhésif, autopolymérisant avec option photopolymérisation
<b>ES Instrucciones de uso</b>	Cemento dental composite autoadhesivo, autopolimerizable con opción de fotopolimerización
<b>PT Instruções de Uso</b>	Composito de cimentação auto-adesivo, autopolimerizável, com opção fotopolimerizável
<b>TR Kullanma Talimatı</b>	Işıkla sertleşme seçeneği kendiliğinden yapışan kendiliğinden sertleşen dental reçine simanı
<b>RU Инструкция по применению</b>	Самоадгезивный самоотверждающий стоматологический композитный цемент с опцией светотверждения

## Try-in of the restoration and isolation

Next, the shade, accuracy of fit and occlusion of the restoration can be checked. Care should be taken when checking the occlusion of fragile and brittle ceramic objects before they are permanently seated, because there is a risk of fracture. It is essential to ensure reliable isolation of the operating field – preferably with OptraGate, cotton rolls, saliva ejector and absorbent pads – when using adhesive cementation with composites.

## Pre-treatment of the abutment surface

According to the manufacturer's directions.

**Note:** Contamination with saliva or blood must be avoided during and after the pre-treatment of the abutment (try-in). If necessary, the abutment should be cleaned again in situ using water and fluoride-free cleaning paste (e.g. with fluoride-free Proxylt).

**A) Abutments made of zirconium oxide and titanium:** Clean or sandblast the abutment surface extraorally as described in section 2.1/2.2

**B) Abutments made lithium disilicate glass-ceramic:** Pre-treat extraorally as described in section 2.3 using e.g. Monobond Etch & Prime

## Sealing of the abutment

- Seal the screw channel of the abutment using e.g. a temporary restorative material (e.g. Telio® CS Inlay). It is recommended to isolate the screw head from the temporary restorative material with cotton pellets.
- Rinse the abutment with water spray.
- Dry the abutment.

## 2 Cleaning and pre-treatment of the restoration

According to the manufacturer's directions.

**Note:** Contamination with saliva or blood must be avoided during and after the pre-treatment of the restoration (try-in).

### 2.1 Restorations made of oxide ceramics - (e.g. zirconium oxide, IPS e.max® ZirCAD)

**IMPORTANT!** In order to achieve a strong bond, do not clean the surfaces with phosphoric acid.

If the restoration has been sandblasted in the laboratory:

- Apply Ivoclean, as described in section 2.4, to clean the restoration.

If the restoration has not been sandblasted in the laboratory:

- Sandblast the inner surface of the restoration (use the blasting parameters stated in the manufacturer's instructions for use of the restorative material, e.g. IPS e.max ZirCAD, max. 1 bar/14.5 psi, Al<sub>2</sub>O<sub>3</sub>).
- Clean the restoration e.g. in an ultrasonic unit for approx. 1 minute.
- Thoroughly rinse the restoration with water spray and dry with water-/oil-free air.

### 2.2 Metal or metal-supported restorations

**IMPORTANT!** In order to create a strong bond, do not clean the metal surfaces with phosphoric acid.

If the restoration has been sandblasted in the laboratory:

- Apply Ivoclean, as described in section 2.4, to clean the restoration.

If the restoration has not been sandblasted in the laboratory:

- Sandblast the inner restoration surface (use the blasting parameters stated in the instructions for use of the restorative material) until an even matt surface has been achieved.
- Clean the restoration e.g. in an ultrasonic unit for approx. 1 minute.
- Thoroughly rinse the restoration with water spray and dry with water-/oil-free air.

**If the restoration comprises precious metal,** apply Monobond® N to the pre-treated surfaces using a brush or microbrush, allow to react for 60 seconds and then disperse with a strong stream of air.

### 2.3 Restorations made of lithium disilicate glass-ceramic (e.g. IPS e.max Press/CAD)

#### A) Procedure using Monobond Etch & Prime

**Regardless of whether the restoration has been pre-treated in the laboratory:**

- Following the try-in, thoroughly rinse the restoration with water spray and dry with water- and oil-free air.
- Apply Monobond Etch & Prime onto the bonding surface using a microbrush and agitate it into the surface for 20 seconds using light pressure. Allow to react for another 40 seconds.
- Then thoroughly rinse off Monobond Etch & Prime with water until the green colour is removed. If residues have remained in the micro-porosities after rinsing, the restoration may be cleaned with water in an ultrasonic cleaner for up to 5 minutes.
- Dry the restoration with a strong stream of water- and oil-free air for approx. 10 seconds.

## B) Procedure using Monobond N

### If the restoration has been pre-treated in the laboratory:

- If the restoration has already been pre-treated in the laboratory, it should be cleaned after the try-in using Ivoclean as described above.
- Subsequently, re-apply Monobond N onto the cleaned surface using a brush or a microbrush, allow to react for 60 seconds and then disperse with a strong stream of air.

### If the restoration has not been pre-treated in the laboratory:

- Etch with 5% hydrofluoric acid (e.g. IPS® Ceramic Etching Gel) for 20 seconds or according to the manufacturer's instructions for use of the restorative material.
- Thoroughly rinse the restoration with water spray and dry with water-/oil-free air.
- Apply Monobond N to the pre-treated surfaces with a brush or microbrush, leave to react for 60 seconds and then disperse with a strong stream of air.

## 2.4 Cleaning of restorations contaminated with blood or saliva

Regardless of any previous conditioning, clean contaminated restoration surfaces extraorally using Ivoclean as follows:

- After try-in, thoroughly rinse the restoration with water spray and dry with oil-free air.
- Cover the entire bonding surface of the restoration with a layer of Ivoclean using a microbrush or brush.
- Allow 20 seconds for Ivoclean to take effect, then thoroughly rinse with water spray and dry with oil-free air.
- Restorations made of lithium disilicate glass-ceramic (e.g. IPS e.max Press/CAD) must be primed again with Monobond Etch & Prime or Monobond N.

## 3 Application of Multilink Speed into the restoration

For each application, place a new automix tip on the double syringe. Extrude Multilink Speed from the automix syringe, apply the desired quantity directly into the restoration and cover the entire bonding surface. As the luting material will cure in the used mixing tip, it may serve as a seal for the contents of the syringe until the next application.

## 4 Placement of the restoration and removal of excess cement

**Note:** As with all composites, Multilink Speed is subject to oxygen inhibition. This means that the surface layer does not polymerize during curing, as it comes in contact with atmospheric oxygen. To avoid this, cover the restoration margins with glycerine gel/air block (e.g. Liquid Strip) immediately after removing the excess material. After complete polymerization, the glycerine gel/air block is rinsed off with water.

### 4.1 Solely self-curing

- Seat the restoration and retain it in place exerting uniform pressure until the excess material has been removed and the cement has completely cured (approx. 3 min.)
- Remove excess material e.g. with a microbrush/brush/foam pellet/dental floss or scaler. Make sure to remove excess material within the processing time, especially in areas that are difficult to reach (proximal or gingival margins, pontics).
- To remove excess material from implant abutments, special implant scalers/curettes should be used to avoid scratching of the implant surface.

### 4.2 Self-curing with additional light-activation to tack-cure excess material or accelerate the curing process (quarter technique, indicated for cases with up to 2 bridge abutments = 3- to 4-unit bridge)

- Seat the restoration and retain it in place exerting uniform pressure until the excess material has been removed and the cement has completely cured.
- Excess cement is light-cured with a polymerization light (approx. 650 mW/cm<sup>2</sup>) from a distance of approx. 0–10 mm for 1 second per quarter surface (mesio-oral, disto-oral, mesio-buccal, disto-buccal).
- Following this procedure, the excess material can be easily removed with a scaler.
- Make sure to remove excess material in time, especially in areas that are difficult to reach (proximal or gingival margins, pontics).
- To remove excess material from implant abutments, special implant scalers/curettes should be used to avoid scratching of the implant surface.
- Subsequently, light-cure all cement lines again for 20 seconds (approx. 1,100 mW/cm<sup>2</sup>). If the light output is lower, curing takes longer.

## 5 Finishing the completed restoration

- Check occlusion and function and adjust as necessary.
- If necessary, rework the cement lines with finishing diamonds.
- Smooth the cement lines using finishing and polishing strips and polish with suitable polishers (e.g. OptraPol®).
- If necessary, finish the restoration margins using suitable polishers (ceramics: e.g. OptraFine®).

## Special instructions for the cementation of endodontic posts:

- For the cementation of endodontic posts, carefully clean the root canal to remove any residual root canal filler material. (Residue of eugenol-based sealers may inhibit the polymerization of the luting composite.)
- Wet the endodontic post, which has been prepared according to the instructions of the post manufacturer, with the mixed Multilink Speed and apply Multilink Speed directly into the root canal using endodontic tips. The material should be applied in a sufficient quantity to provide excess cement.
- Insert the endodontic post.
- Remove excess cement from the occlusal preparation surface. Then, light-cure Multilink Speed for 20 seconds.
- Pre-treat the occlusal preparation surface with an adhesive (e.g. Tetric® N-Bond Universal) according to the respective Instructions for Use.
- Apply the core build-up material (e.g. MultiCore®) directly onto the occlusal preparation surface and cure according to the manufacturer's instructions (light-curing).

## Special notes

Multilink Speed should have room temperature when it is applied. At refrigerator temperature, the material is difficult to extrude and its processing and curing times will be longer.

Do not clean the Multilink Speed automix syringes with agents that have an oxidative effect.

# CE 0123

Rx ONLY

**Date information prepared:**  
2019-06-03/Rev. 0  
702853/7spr

**Manufacturer**

Ivoclar Vivadent AG  
FL-9494 Schaan/Liechtenstein

## Warning

- Avoid contact of uncured Multilink Speed with the skin/mucous membrane or eyes.
- Uncured Multilink Speed may cause slight irritation and may lead to a sensitization against methacrylates.
- Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.

## Storage

- Store Multilink Speed at 2–28 °C/36–82 °F.
- Do not use Multilink Speed after the indicated expiry date.
- Do not disinfect syringes with oxidizing disinfectants.
- Store Multilink Speed automix syringe with the mixing tip attached after use.
- Expiry date: see note on syringes and packages.

## Keep out of the reach of children!

### For use in dentistry only.

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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