# Centric Tray

Verarbeitungsanleitung / Instructions for Use / Mode d'emploi  
Istruzioni d'uso / Instrucciones de uso / Instruções de Uso

CE

ivoclar  
vivadent®  
technical

<b>Gebrauchsinformation</b>	<b>Seite 2</b>
<b>Instructions for Use</b>	<b>Page 12</b>
<b>Mode d'emploi</b>	<b>Page 22</b>
<b>Istruzioni d'uso</b>	<b>Pagina 32</b>
<b>Instrucciones de uso</b>	<b>Página 42</b>
<b>Instruções de Uso</b>	<b>Page 52</b>

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# Centric Tray with Stratos 200 Articulator

# Centric Tray

The success of a dental prosthetic restoration largely depends on the accurate determination of centric relation (position of the lower jaw in relation to the upper jaw) in the patient. Since the upper and lower jaw in centric relation constitute a closed three-dimensional unit in the mouth, individual impressions of the jaws cannot capture the totality of the relationship. Ivoclar Vivadent offers the Centric Tray to determine the occlusal position. This tray permits the technician to mount existing models in the proper relation – a prerequisite for accurate placement of intraoral registration pins. Additionally, the keyed handle on the Ivoclar Vivadent Centric Tray enables the cranial/articulation oriented transfer with the UTS Transferbow in the articulator.

## Indication


- a) Edentulous patients
- b) Partially edentulous patients

## Advantages

- Impressions / bite in one appointment
- Easy, accurate location of occlusal position
- Easy assembly of the models from separate impressions
- Universally applicable to most jaw sizes
- Supports can be placed individually
- Lip support
- UTS Transferbow can be attached
- Physiological alignment of intraoral centre pin recordings (e.g. Gnathometer "M")
- Autoclavable and reusable (It is advisable to autoclave the tray not more than 10 times. If trays are used which have been autoclaved more than 10 times, the stability of shape required for accurate dimension recording is no longer ensured.)

## The Centric Tray and its components


- 1. Centric Tray
- 2. UTS Transferbow connection
- 3. Lip support
- 4. Retaining track
- 5. Regular support
- 6. Supporting wings

# Clinical Application

## a) Edentulous cases

- 1) Take separate impressions, as normal.
- 2) Slip the support wings (6) distally onto the retaining track (4). Try-in the Ivoclar Vivadent Centric Tray, and adjust by trimming as necessary.


- 3) To record the occlusal height, first make marks on the patient's nose and chin, using a suitable marker. Then proceed by selecting one of the options below.
  - Replace the patient's dentures and measure the present vertical dimension
  - Position the empty Centric Tray in the mouth. Ensure that the lips are moist and seated on the lip support (3), and that they gently support the tray. Leave the tray in-situ for 1–2 minutes, then measure and record the vertical dimension. Remove tray from mouth.

- 4) Mix silicone impression material, alginate or similar material and load both sides of the tray.

Recommended trays: Accutrax System 1 for impression taking of jaws; distribution in Europe by Candulor (Switzerland).

Tray with alginate


Tray without alginate

5) Insert the loaded tray into the mouth sideways, rotating into position over the lower ridge. Place the moist lips on the lip support (3), whilst slowly closing the mouth until the marked occlusal height is reached. Allow the silicone, alginate or other material used to set completely. **Check the vertical measure** and remove the tray.

6) The UTS Transferbow conveys patient-related information to the articulator. Align the Transferbow to the correct coordinates and attach it to the Ivoclar Vivadent Centric Tray at the UTS connection point (2). Fix the registration joint and remove the Transferbow.

7) The occlusal record, with the UTS Transferbow attached, is now handed to the technician, who will set-up the teeth according to the patient's coordinates.

#### Notes:

- Keep a record of the measured occlusal height.
- Use a viscous silicone, alginate or similar material.
- Use sufficient impression material.


## b) Partially edentulous cases

1) Take separate impressions, as normal.

We recommend using the Accutray System 2 for impression taking of the jaws; distribution in Europe by Candulor (Switzerland).

Tray with alginate

Tray without alginate

2) Check the closed occlusal position. If, due to the loss of contact between the opposing antagonists this is not possible, proceed as in Edentulous Cases (3).

3) Slip the regular supports (5) distally onto the retaining track (4), placing them into the gaps. Try-in the Ivoclar Vivadent Centric Tray and adjust by trimming as necessary.


4) Mix silicone impression material and load both sides of the tray.

5) Insert the loaded tray into the mouth sideways, rotating into position over the lower ridge. Place the moist lips on the lip support (3), whilst slowly closing the mouth until the marked occlusal height is reached. Allow the silicone, alginate or other material used to set completely. **Check the vertical measure** and remove the tray.

6) The UTS Transferbow conveys patient-related information to the articulator. Align the Transferbow to the correct coordinates and attach it to the Ivoclar Vivadent Centric Tray at the UTS connection point (2). Fix the registration joint and remove the Transferbow.

7) The occlusal record, with the UTS Transferbow attached, is now handed to the technician, who will set-up the teeth according to the patient's coordinates.


# Technical Instructions

1) Trim back the vestibular edges of the silicone impression. Try to identify three reference points to the models.

2) Seat the models in the silicone impression and check the fit.

3) Mount the models in the articulator, positioning them in relation to the occlusal record.

- 4) If a UTS Transferbow has been used, attach it to the articulator, applying the philosophy relative to that determined by the Transferbow record.

- 5) Secure the Ivoclar Vivadent Centric Tray in place. Place and fix the maxillary model. Remove lower jig and locate and fix the mandibular model.
- 6) Fabricate individual trays, centre pin recording devices, wax occlusal rims or other prescription requests.

#### **Note**

- Check the seating of the models.
- Consider the articulator philosophy.

# Delivery form

- 1 Centric Tray
- 3 Regular supports
- 3 Supporting wings
- 1 Set of Instructions for Use

## **Ivoclar Vivadent Centric Tray – a complementary product of BPS®**

The Ivoclar Vivadent Centric Tray is part of the future-oriented dental prosthetic system, which includes individualized dental lab marketing activities. You too can benefit from the advantages of the Biofunctional Prosthetic System (BPS) in your day-to-day work.

This tray has been developed solely for use in dentistry. Assembly and operation must be undertaken strictly in accordance with the Instructions for Use. Liability for damages will not be accepted as a result of misuse or any failure to observe the Instructions. The user is solely responsible for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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## **Ivoclar Vivadent Limited**
## **Ivoclar Vivadent, Inc.**